EN FR DE IT ES TR

CE

0123

MD

Gel Desensitizer

Desensitizer for Hypersensitive Teeth

EN

Instructions for Use

IMPORANT : READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE. KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Desensitizer for Hypersensitive Teeth

Gel Desensitizer

FOR DENTIST USE ONLY

1. Product overview

Gel Desensitizer is a single-syringe dental desensitizer with fluoride. It is intended to be applied to tooth surfaces for treatment of dental hypersensitivity. Hypersensitivity can be caused by various reasons such as toothbrush abrasion, gingival recession, periodontal disease and acid erosion. Gel Desensitizer reacts chemically with calcium in teeth and forms precipitates which can seal not only dentinal tubules, but also microcracks before and/or after whitening. Gel Desensitizer suppresses hypersensitivity.

Intended purpose:

Gel Desensitizer is intended to treat dental hypersensitivity.

Indication:

Dental hypersensitivity on exposed dentin and/or microcracks in the enamel.

Contraindication:

Gel Desensitizer should not be used in patients who have a history of hypersensitivity reaction (e.g. allergy, inflammation, ulcer, and rash) to this product.

Possible side effect:

Gel Desensitizer or one of its components may cause hypersensitive reactions (e.g. allergy, inflammation, ulcer, and rash).

2. Components

Gel Desensitizer

① Gel Desensitizer ----- 3mL ----- Class IIa

CE

0123

CE

② Plastic Needle Tip ----- 10pcs ----- Class I

■ Attached document: Instructions for Use

3. Composition

Gel Desensitizer: Copolymer with sulfonic acid group, Oxalic acid, Water, Sodium fluoride, pH buffer (Potassium salt), Thickener and the others.

Component	Physical Form
Gel Desensitizer	Pale yellow to milky-white viscous liquid

4. Precautions

This product contains substances that may cause allergic reactions. Read the following points.

1) Safety precautions

- ① If accidentally swallowed during use, immediately have the patient gargle with water, consult a physician if necessary.
- ② Discontinue the use immediately if signs of irritation (rash, etc.) appear. See a physician.
- ③ This product is acidic (approx. pH 2-3). Avoid contact with wound or inflammation on the oral mucosa or soft tissue. If this product contacts wound or inflammation on the oral mucosa or soft tissue, immediately rinse the contacted area with copious amounts of water. In addition, avoid contact with eyes. In case of eye contact, flush immediately with copious amounts of water and consult an ophthalmologist if necessary. Residual Gel Desensitizer on the lip or the corner of the mouth may cause hypersensitivity reactions. If a patient shows hypersensitivity reactions, flush the applied area with copious amounts of water and have the patient seek for a medical care.

2) Common precautions

- ① When using this product, the operator should wear protective gloves and eyewear.
- ② If an operator has a history of hypersensitivity to this product, use protective gears to avoid direct contact to it. If the hypersensitivity reactions appear, see a physician.
- ③ Be cautions on a risk of causing gingival inflammation by rubbing the gingiva excessively with cotton pledget and rubber cup etc.
- ④ Do not mix components with other materials.
- ⑤ Use this product only for the application recommended in these instructions.
- ⑥ When using this product in combination with other dental materials or instruments, the usage and precautions described in the instructions attached to each product should be strictly observed.
- ⑦ Use common agents to clean and disinfect the syringe in accordance with the commonly used hygienic regulations.

3) Handling precautions

- ① Desensitizing efficacy may be reduced if the procedures in this instruction is not followed.
- ② Plastic Needle Tip is one-time use only. Dispose it after each patient use to avoid cross-contamination.
- ③ Tightly close syringe with syringe cap immediately after use to prevent the material from drying out. Plastic Needle Tip will not serve as a syringe cap.
- ④ Gel Desensitizer should be used only with Plastic Needle Tip provided.
- ⑤ If you push the plunger strongly, the gel might splatter around. Please push the plunger gently.
- ⑥ Do not use Gel Desensitizer after the root planing if reattachment is expected.
- ⑦ Desensitizing efficacy may be reduced if a patient has bruxism or a tooth with a large non-carious cervical lesion (wedge-shaped defect).
- ⑧ Once a tooth is sealed with Gel Desensitizer, the bond strength of subsequently applied adhesives will be compromised.

- If a surface treated with Gel Desensitizer must be bonded, expose the fresh surface by sanding or etching. Application of Super-Bond C&B Green Activator or Dentin Etchant Gel (manufactured by Sun Medical Co., Ltd.), rubbed for 10 seconds with a saturated pledget before bonding with Super-Bond (manufactured by Sun Medical Co., Ltd.), will attain good bond strength.
- ③ This product contains fluoride. Avoid contact with the glass product as it can react/etch.
- ④ Material's taste may change slightly depends on the storage environment and duration. However, it does not affect the performance.

5. Storage

1) Storage conditions

- ① Store Gel Desensitizer in dark at room temperature (1°C - 30°C). Avoid exposure to direct sunlight and extreme temperatures (above 30°C or below 1°C).
- ② Do not use this product after the expiration date. See outer package for expiration date.

2) Product lifetime

Product lifetime can be shortened depending on the circumstances of use. Read instructions of relevant products carefully before use.

6. Disposal

Dispose of contents, containers, Plastic Needle Tip, and packages in accordance with local regulations.

As in any dental treatments, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditons for use.

7. How to use

① **Surface preparation**

Tooth surface preparation: Remove plaque and clean the tooth surface with a cotton pledget and floss thoroughly.

② **Application of Gel Desensitizer**

- Remove the syringe cap by twisting it counter-clockwise. Save the syringe cap.
- Install the Plastic Needle Tip to the syringe properly and make sure the Plastic Needle Tip and the syringe are tightly connected.
- Apply a thin layer to the lesion or the tooth surface around it and leave it for 30 seconds. Recommended quantity per tooth is the size of a grain of rice (φ:3mm).
- Dispose Plastic Needle Tip after use and seal the container immediately with the syringe cap. Do not store the syringe with Plastic Needle Tip in place. Store it only with the syringe cap.

③ **Rinse**

Have the patient rinse the mouth with water thoroughly.

If the desensitizing is not satisfactory with the first application, repeat the treatment at successive appointments to increase desensitization.

* To obtain higher desensitizing results or for the case establishing the clean surface is difficult due to the severe sensitivity, apply Gel Desensitizer with **rubbing/pumping motion for more than 5 seconds** using cotton pledget or rubber cup by following the instruction below.

- Selection and preparation of applicators
 - Select a cotton pledget according to the size of the affected area or attach the rubber cup etc. to the contra-angle.
- Rubbing/pumping motion
 - Dispense the desensitizer onto the lesion and/or surrounding tooth surface (standard quantity per tooth is the size of a grain of rice). Rub and pump for more than 5 seconds using the cotton pledget. If using the rubber cup, pick up the desensitizer (about the size of a grain of rice per tooth) in the cup, apply by rubbing and pumping for more than 5 seconds at low speed (around 1,000rpm) so that Gel Desensitizer does not scatter. The use of cotton pledget and rubber cup etc. is preferable for the lesions close to gingiva or interproximal areas.

When any serious incident has occurred in relation to the device, please report to the following manufacturer and the competent authority in your country.

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