# Gel Desensitizer

## Desensitizer for Hypersensitive Teeth

EN

### INSTRUCTIONS

**IMPORTANT :**  
 READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.  
 KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Desensitizer for Hypersensitive Teeth

## Gel Desensitizer

FOR DENTIST USE ONLY

### 1. What is Gel Desensitizer?

Gel Desensitizer is a syringable dentin tubule sealant. Copolymer with a sulfonic acid group and oxalic acid in Gel Desensitizer react chemically with calcium in dentin / enamel to form a protective polymer coating containing calcium oxalate crystals. This coating also contains fluoride and potassium salt, it tightly seals not only the tooth surface and dentinal tubules but also the internal micro-cracks reaching to dentin. Sealing the dentinal tubules can suppress hypersensitivity symptom by reducing the fluid movement in the tubules, which is one of the etiologies of the hypersensitivity.

Gel Desensitizer is indicated for topical treatment of dental hypersensitivity. The use of Gel Desensitizer will NOT relieve the discomfort of pulpitis.

### 2. Contents

Gel Desensitizer

- ① Gel Desensitizer ----- 3mL
- ② Plastic Needle Tip ----- 10 (with 1 grey Needle Cap)
- ③ Instructions

### 3. Precautions : Read all instructions thoroughly before use.

#### 3.1 Safety

##### 1) Precautions for use

- ① Avoid mixing with other products.
- ② If accidentally swallowed during use, immediately have the patient gargle with a water, consult a physician if necessary.
- ③ Plastic needle is one-time use only.

##### 2) Important Precautions

- ① Discontinue the use immediately if signs of irritation (rash, etc.) appear. See a physician.
- ② If an operator has a history of hypersensitivity to this material, use protective gears to avoid direct contact to it. If the hypersensitivity reactions appear, see a physician.
- ③ This product is acidic. Avoid contact with oral mucosa, skin or eyes. Contaminated skin should be thoroughly rinsed with running water. Materials left on the tissue may cause, for example, swelling of the lips or white sore at the corner of the mouth. In case of the eye contact, flush immediately with copious amounts of water and consult an ophthalmologist if necessary.
- ④ Be cautious on a risk of causing gingival inflammation by rubbing the gingiva excessively with cotton pledget and rubber cup etc.
- ⑤ Make the best clinical judgement on selecting the treatments and materials, as each patient and case has its unique requirements.

##### 3) Warning of adverse events

The use of this material may produce allergic skin reactions in certain sensitive individuals.

#### 3.2 To get the best results with Gel Desensitizer

##### 1) Storage conditions

- ① Store Gel Desensitizer in dark at room temperature. Avoid exposure to direct sunlight and extreme temperatures (above 30°C or below 1°C).
- ② Store securely with access only to the dental professionals.

##### 2) Precautions about how to use

- ① Do not use Gel Desensitizer after the root planning if reattachment is expected.
- ② Desensitizing efficacy may be reduced if a patient has bruxism or a tooth with a large non-carious cervical lesion (wedge-shaped defect).
- ③ Once a tooth is sealed with Gel Desensitizer, the bond strength of subsequently applied adhesives will be compromised.  
 To seal a surface that will be bonded later, use Hybrid Coat (or Hybrid Bond) instead of Gel Desensitizer.  
 If a surface treated with Gel Desensitizer must be bonded, expose the fresh surface by sanding or etching. Application of Green Activator or Dentin Etchant Gel (sold separately), rubbed for 10 seconds with a saturated

pledget before bonding with Super-Bond C&B, will attain good bond strengths.

- ④ This material contains fluoride. Avoid contact with the glass product as it can react/etch.
- ⑤ Recap the syringe immediately after each use. Dried out material can clog the syringe if left uncapped.
- ⑥ Material's taste may change slightly depends on the storage environment and duration. However, it does not alter the performance.

### 4. How to use Gel Desensitizer

#### ① Surface preparation

Tooth surface preparation: Remove plaque and clean the tooth surface with a cotton pledget and floss thoroughly.

#### ② Application of Gel Desensitizer

- a) Remove the syringe cap and attach a disposable plastic needle to the syringe.
- b) Apply a thin layer to the lesion or the tooth surface around it and leave it for 30 sec. Recommended quantity per tooth is the size of a grain of rice (Φ:3mm).

#### ③ Rinse

Have the patient rinse the mouth with water thoroughly.

***If the desensitizing is not satisfactory with the first application, repeat the treatment at successive appointments to increase desensitization.***

\* To obtain higher desensitizing results or for the case establishing the clean surface is difficult due to the severe sensitivity, apply Gel Desensitizer with **rubbing/pumping motion for more than 5 seconds** using cotton pledget or rubber cup by following the instruction below.

##### 1) Selection and preparation of applicators

Select a cotton pledget according to the size of the affected area or attach the rubber cup etc. to the contra-angle.

##### 2) Rubbing/pumping motion

Dispense the desensitizer onto the lesion and/or surrounding tooth surface (standard quantity per tooth is the size of a grain of rice). Rub and pump for more than 5 seconds using the cotton pledget. If using the rubber cup, pick up the desensitizer (about the size of a grain of rice per tooth) in the cup, apply by rubbing and pumping for more than 5 seconds at ~1,000 rpm so that Gel Desensitizer does not scatter. The use of cotton pledget and rubber cup etc. is preferable for the lesions close to gingiva or interproximal areas.

# Gel Desensitizer

## Desensitizer for Hypersensitive Teeth

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#### Gel Desensitizer

- ①Gel Desensitizer ----- 3mL
- ②Plastic Needle Tip ----- 10 (with 1 grey Needle Cap)
- ③Instructions

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SUN MEDICAL CO.,LTD.