BONDFILL DUAL

EN

INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Dental Adhesive Resin for Filling
 Polymer-based Pit and Fissure Sealant, Class 1

Bondfill DUAL

FOR DENTIST USE ONLY

1. What is Bondfill DUAL?

Bondfill DUAL is a dual-cure, self-adhesive resin for dental restorations. It contains 4-META/MMA-TBB. Though it requires no bonding agent, Bondfill DUAL forms strong bonds to enamel, dentin, and restorative materials without creating gap and discoloration. It is also designed to absorb intense, complicated complex external stresses and has low wear rate. The Powder contains PMMA and TMPT reactive organic fillers which creates a tough restoration that is hard yet still very resilient. The inorganic filler loading is 1.2 vol% (as powder) with a particle size range of 5 to 50 nm.

Typical Applications:

Bondfill DUAL is indicated wherever resilience is required, such as:

- Restoration of tooth wear, such as attrition, abrasion, abfraction, wedge-shaped defects
- Secondary caries at a restoration margin
- Repair of fractured prostheses and damaged restorations
- Sealing incipient pit and fissure caries

Note: Bondfill DUAL is not indicated to use for the restoration of occlusal surface.

2. Contents of the kit

2.1 Bondfill DUAL Kit:

- ①Dual Base ----- 1 (8mL)
- ②Catalyst V ----- 1 (0.7mL)
- ③Powder Cervical (approx. Shade A4) ----- 1 (3g)
- ④Powder Medium (approx. Shade A3.5) ----- 1 (3g)
- ⑤Teeth Primer ----- 1 (3mL)
- ⑥Sponge (L•S) ----- 1
- ⑦Dispensing Stand ----- 1
- ⑧Dispensing Cups ----- 20
- ⑨Brush Handle (Curved) ----- 1
- ⑩Brush Tips (Brush-dip L) Label color: Pink----- 10
- ⑪Brush Tips (Brush-dip LL) Label color: Purple ----- 10

Attached documents : Instructions / Pictorial Instruction Card

2.2 Refills:

- Dual Base ----- 8mL
- Catalyst V -----0.7mL
- Powder Light (approx. Shade A2) ----- 3g
- Powder Medium (approx. Shade A3.5) ----- 3g
- Powder Cervical (approx. Shade A4) ----- 3g
- Powder Opacious (approx. Shade A3 Opaque)----- 3g
- Teeth Primer----- 3mL
- Sponge (L•S)----- 1
- Dispensing Stand ----- 1
- Dispensing Cups ----- 40
- Brush Handle (Straight) ----- 1
- Brush Handle (Curved) ----- 1
- Brush Tips (Brush-dip S) Label color: Green ----- 10
- Brush Tips (Brush-dip L) Label color: Pink ----- 10
- Brush Tips (Brush-dip LL) Label color: Purple----- 10

3. Precautions : Read all instructions thoroughly before use.

3.1 Safety

- Sensitivity**
 Bondfill DUAL should not be used by clinicians or on patients who are sensitive to methacrylates or acetone. If signs of irritation such as redness appear, stop using Bondfill DUAL immediately and consult a physician.
- Avoid contact**
 Avoid contacting soft tissues, skin or eyes by Bondfill DUAL components. The use of rubber dam is recommended. Contaminated skin or mucosa should be wiped immediately with alcohol and then rinsed thoroughly with copious amounts of running water. If Bondfill DUAL components enter the eye, rinse immediately and thoroughly with running water and have the patient examined by an ophthalmologist. Clinicians should use dental gloves.
- Flammability**
 Catalyst V, Dual Base, and Teeth Primer are flammable. Do not store where they may be exposed to open flame.
- Applications**
 Use Bondfill DUAL only for the applications recommended in these instructions.

3.2 Precautions on Catalyst V

Catalyst V reacts with air and water to generate heat and lose activity.

Please abide by the following.

- Storage conditions**
 Avoid high temperature, high humidity and direct sunlight. The repeated temperature changes may cause the syringe to aspirate air, thereby shortening the Catalyst V's shelf-life.
 *After prolonged storage, the first drop of the Catalyst V may be inactive due to exposure to air, though the rest of the material remains active. Express and dispose of a drop before making the mix.
 *The syringe is made of glass, therefore it must be handled with care to prevent shock, dropping, and other physical damage.
- Cap closure**
 The cap simply slides on and off. Recap the syringe immediately after each use. Air (oxygen and humidity) deactivates the Catalyst V. Do not leave the cap off during the bonding procedure.
- Expressing the Catalyst**
 If the Catalyst V does not come out of the syringe because the tip is blocked or the screw is tight, do not try too hard to turn it. The contents may splash as the syringe breaks.
- After use**
 Unscrew the syringe screw two turns counter-clockwise after each use to relieve pressure on the Catalyst. (Pressure buildup can cause leakage of the Catalyst V or even crack the syringe.)
- Clean spilled Catalyst V immediately with wet towels**
 The Catalyst V reacts with oxygen. If absorbed by a flammable material, it may raise the temperature enough to cause smoldering. If the Catalyst V is spilled, wipe it up immediately with a WET (not dry) disposable towel. Then rinse the towel to deactivate the Catalyst V thoroughly in running water.
- Cleaning syringe**
 Wipe the tip of the syringe with a dry gauze after each use to prevent residue buildup. Then rinse the gauze with water to kill any remaining activity. Buildup of the residue may prevent the cap from seating properly.

3.3 Storage

Please take the following precautions to maintain the quality.

- Storage conditions**
 Store Catalyst V, Dual Base and Teeth Primer in a cool (1°C - 30°C) and dark location. High temperature, high humidity and direct sunlight will shorten their shelf-life.
- Volatility**
 Teeth Primer and Dual Base are highly volatile. Recap the bottle immediately.
- Contamination**
 To avoid contamination, do not mix the components with other materials. Firmly recap the containers immediately after use. Discard any excess of mixed materials after use. Once dispensed, do not return materials to the containers. Both Dispensing Cups and Brush Tips are for one-time use. Dispose of them after use.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

4. How to use Bondfill DUAL

a. Surface preparation

- Cleaning**
 Remove contaminants and stains using a polishing bur or oil-free, fluoride-free pumice to expose fresh tooth surface. Rinse thoroughly and dry. A rubber dam is highly recommended to reduce the chance of contamination by saliva, breath humidity or blood.
- Priming**
 Dispense Teeth Primer into one of the small recesses on the Dispensing Stand. Using a pledget, keep the surface to be bonded saturated with a copious amount of Teeth Primer for approximately 20 seconds. Care should be taken not to contact the gingiva with the pledget. Then air dry. Rinsing is not needed.

Standard application time

Surface	Teeth Primer	Enamel Etchant Gel	Dentin Etchant Gel
Enamel	Approx. 20 sec	30 sec	30-60 sec
Dentin		—	5-10 sec
After application	Dry without rinsing	Rinse and dry	

*If the surface is cut enamel, use either Teeth Primer or Enamel Etchant Gel. If the surface is uncut enamel, however, Enamel Etchant Gel is recommended. If the surface is dentin, use either Teeth Primer or Dentin Etchant Gel.
 *Both Enamel Etchant Gel and Dentin Etchant Gel are sold separately by Sun Medical Co., Ltd. Unlike the Teeth Primer (which should not be rinsed from the tooth surface), Enamel Etchant Gel and Dentin Etchant Gel must be rinsed off after the proper etching period (see table above).

*If a precious metal or ceramic restoration is present on the tooth surface to be treated, prime the metal or ceramic surface after conditioning the tooth surface with Teeth Primer.

3) Metal Surface

Rinse the metal surface thoroughly and dry. Applying the Enamel Etchant Gel then rinse after can also help in removing surface contaminants. Air-abrade the metal surface with 50-micron aluminum oxide. If the surface is precious metal, apply a metal primer such as V-PRIMER on the metal surface to be bonded in order to increase bond strength and durability. V-PRIMER is sold separately by Sun Medical Co., Ltd. For application details refer to the instructions for the metal primer.

4) Ceramic (porcelain and zirconia) surface

Depending on the case at hand, clean the porcelain/zirconia surface with abrasive wheels and points or a prophy cup with fluoride-free, oil-free pumice.

Rinse thoroughly with water and dry. Applying the Enamel Etchant Gel helps in removing surface contaminants. Then, apply a porcelain primer such as Porcelain Liner M, Universal Ceramic Primer or Zirconia Liner Premium on the porcelain/zirconia surface to be bonded in order to increase bond strength and durability. These ceramic primers are sold separately by Sun Medical Co., Ltd. Refer to the instruction of each product for the application details.

b. Preparation of Bondfill DUAL

Place two Dispensing Cups into two of the three holes in the Dispensing Stand. Dispense 2-3 drops of the Dual Base into one of the Dispensing Cups. Hold the Catalyst V syringe vertical, and turn the screw to dispense a drop into the Dual Base. Stir gently with a brush. This mixture is called "activated liquid." The activated liquid gradually loses its potency. Use it within 3 minutes after the activation.

Activated Liquid

Intended Use	Dual Base	Catalyst V	Operartion time
Restoration and Build-up	2-3 drops	1 drop	3 minutes after mixing
Pit and Fissure Sealant	2 drops		

c. Restoration/Build-up/Pit and Fissure Sealant

- Dispense the appropriate amount of Powder into the other Dispensing Cup.
- Wet the surface to be bonded with the activated liquid.
- Soak the brush with ample activated liquid, then dip the brush into the Powder to form a ball of powder on the tip. Brush the powder ball onto the pre-wet surface to be bonded.
- Repeat the procedure until the restoration or the build-up is completed. When you repeat the procedure, clean the brush tip with gauze between dips.
- Use LED/Halogen light irradiator for 20 seconds or Plasma arc light irradiator for 6 seconds as necessary. The polymerization process can be accelerated.

d. Setting and finishing

For better results, Sun Medical recommends:

- If curing with light irradiator (LED/Halogen: for 20 seconds) or (Plasma arc : for 6 seconds), wait approximately 7 minutes after placing the last powder ball before finishing and polishing on the same day.
- In the case of no irradiation : wait approximately 10 minutes before finishing and polishing.
- If complete setting is desirable, wait to polish until the next day.

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Dental Adhesive Resin for Filling
Polymer-based Pit and Fissure Sealant, Class 1

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2.1 Bondfill DUAL Kit:

- ①Dual Base ----- 1 (8mL)
- ②Catalyst V ----- 1 (0.7mL)
- ③Powder Cervical (approx. Shade A4) ----- 1 (3g)
- ④Powder Medium (approx. Shade A3.5)----- 1 (3g)
- ⑤Teeth Primer ----- 1 (3mL)
- ⑥Sponge (L·S) ----- 1
- ⑦Dispensing Stand ----- 1
- ⑧Dispensing Cups ----- 20
- ⑨Brush Handle (Curved) ----- 1
- ⑩Brush Tips (Brush-dip L) Label color: Pink ----- 10
- ⑪Brush Tips (Brush-dip LL) Label color: Purple ----- 10

Attached documents : Instructions / Pictorial Instruction Card

2.2 Refills:

- Dual Base ----- 8mL
- Catalyst V ----- 0.7mL
- Powder Light (approx. Shade A2) ----- 3g
- Powder Medium (approx. Shade A3.5)----- 3g
- Powder Cervical (approx. Shade A4) ----- 3g
- Powder Opacious (approx. Shade A3 Opaque)----- 3g
- Teeth Primer ----- 3mL
- Sponge (L·S) ----- 1
- Dispensing Stand----- 1
- Dispensing Cups ----- 40
- Brush Handle (Straight) ----- 1
- Brush Handle (Curved) ----- 1
- Brush Tips (Brush-dip S) Label color: Green ----- 10
- Brush Tips (Brush-dip L) Label color: Pink ----- 10
- Brush Tips (Brush-dip LL) Label color: Purple ----- 10

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