## INSTRUCTIONS

IMPORTANT :  
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.  
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Dental Adhesive Resin Cement

# Chemiace II

FOR DENTIST USE ONLY

# Chemiace II

## 1. What is Chemiace II?

Chemiace II is a dual-cure resin cement for onlays, inlays, crowns and bridges. It contains 4-META\*, a high performance adhesive monomer that creates excellent bonds to dentin and enamel by virtue of formation of a hybrid layer. Handling is similar to glass ionomers, while chemically Chemiace II is a resin cement and thus water resistant. Removal of excess resin is quick and easy. The cement is radiopaque. The pretreatment with V-PRIMER for precious metals and Porcelain Liner M for porcelain will further enhance the bond strength and its durability.

\* 4-methacryloxyethyl trimellitate anhydride

## 2. Contents of the kit:

- ① Powder ----- 7g
- ② Liquid ----- 6.5mL
- ③ Green Activator ----- 5mL
- ④ Spatula ----- 1
- ⑤ Measuring Spoon ----- 1
- ⑥ Mixing Pad ----- 1

## 3. Precautions : Read all instructions thoroughly before use.

### 3.1 Safety

#### ① Sensitivity

Chemiace II should not be used by clinicians or on patients who are sensitive to methacrylates or resin cements. If signs of irritation such as redness appear, stop using Chemiace II immediately and consult a physician.

#### ② Avoid contact

Avoid contacting soft tissues, skin or eyes by Chemiace II components. The use of rubber dam is recommended. Contaminated skin or mucosa should be wiped immediately with alcohol and then rinsed thoroughly with copious amounts of running water. If Chemiace II components enter the eye, rinse immediately and thoroughly with running water and have the patient examined by an ophthalmologist. Clinicians should use dental gloves.

#### ③ Pulp protection

If the preparation approaches the pulp, apply a protective liner or base.

#### ④ Flammability

The Liquid is combustible. Keep away from open flame in use and storage.

#### ⑤ Give care to acidity

The Green Activator is acid. Do not make patients swallow it during use.

#### ⑥ Applications

Only use Chemiace II for the applications recommended in these instructions.

### 3.2 To get the best results with Chemiace II

#### ① Storage conditions

Store Chemiace II in a cool, dark place. **LIQUID should be stored IN A REFRIGERATOR.** High temperature and direct sunlight will shorten their shelf-life.

#### ② Avoid Contamination

To avoid contamination, do not mix the components with other materials. Firmly recap the containers immediately after use. After the cementing procedure, dispose of any excess of mixed materials. If the mixture is spread into a thin film, it will remain usable for about 10 minutes. Once dispensed do not return materials to the container.

**As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.**

## 4. How to use Chemiace II

### ① Surface preparation

#### 1) Tooth Surface

Clean cavity/preparation, wash with water and dry with oil-free air. (A rubber dam is highly recommended to reduce the chance of contamination by saliva, breath humidity or blood.) Then apply the Green Activator to the tooth surface. After the appropriate treatment time (see table), rinse well with water. Dry the entire surface with air syringe.

	Treatment time
Enamel	30-60 seconds
Dentin	5-10 seconds

\*If both dentin and enamel are involved, apply the Green Activator first to the enamel, and after 30 seconds to the dentin. After an additional 5-10 seconds, rinse both surfaces with water and air dry.

### 2) Restoration Surface

To the extent practically possible, the surface of a restoration to be bonded should be abraded using a sand-blast or roughened by a round bur. Rinse thoroughly and dry. Then, depending on the restorative material, the roughened and cleaned surface may be treated with a surface modifier, such as Porcelain Liner M for porcelain or V-PRIMER for precious metals, to get higher bond strength and durability.

### ② Mixing and Application of the cement

- For a standard mix, using the Measuring Spoon dispense one scoop (small cup) of the Powder onto the mixing pad. For the best mixing ratio, it is recommended the cup is overfilled with powder. The powder is then leveled by pressing the spoon firmly against the leveling bar situated at the top of the jar while pulling in a horizontal direction. (The Powder weight scooped with the small cup should be approx. 0.03g.)
- Then dispense one drop of the Liquid onto the pad. (For accurate liquid dispensing, hold the bottle vertically to reduce air bubbles and squeeze to express a drop. Further, the bottle nozzle should be wiped after use to prevent residue buildup.) Both Powder and Liquid containers should be recapped immediately after use to prevent possible contamination.
- The dispensed Powder and Liquid should be mixed well using the Spatula. The flow and the setting time of mixed Chemiace II can be adjusted by altering the Powder/Liquid ratio.

	Powder / Liquid ratio	Setting time	Flow (mm)	Bond strength (MPa)
Standard ratio	1 scoop (small cup) / 1 drop (standard quantity for an inlay)	4.5 min	55	12
	1 scoop (big cup) / 3 drops (standard quantity for a crown)			
High	1 scoop (big cup) / 2 drops	4 min	35	12
Low	2 scoops (small cup) / 3 drops	5.5 min	60	12

\*3 scoops of small cup = 1 scoop of big cup

- Apply the mixture to the surfaces of both the prepared tooth and the restoration.
- Note:** The mixture should be used within 1 minute of mixing or its setting becomes slow. If the mixture is spread into a thin film so that its polymerization is oxygen-inhibited, use can be delayed for up to 10 minutes.

### ③ Cementing, setting and finishing

- Insert the restoration immediately after the application of the cement. After confirming that it is completely seated, hold it under pressure in position until the cement sets.
- To remove the excess cement, irradiate the excess with visible light for 10 seconds and immediately remove it. Or, alternatively without light, remove the excess on the halfway to setting, specifically 1.5 - 2 minutes after cementing.
- The final setting time with the standard Powder/Liquid ratio is 4.5 minutes in the mouth (when used immediately after mixture), thereafter polish the cement line or adjust occlusion according to normal practices.

#### Option: Light-setting

If the restoration is translucent like porcelain or cured composite resin, light-setting may be opted. In this case, irradiate visible light first for 10 seconds, remove the excess cement, and then irradiate the light another 30 seconds for the final setting.

## 5. Data

### Bond strength to various materials

Material to be bonded	Pre-treatment	Tensile Bond Strength (MPa)	
		TC 5,000	TC 10,000
Bovine dentin	Green Activator 10sec	12	
Bovine enamel	Green Activator 30sec	10	
Gold alloy (Type IV)	Sandblasting	12	3
		14	7
Pd/Ag/Au alloy	Sandblasting + V-PRIMER	13	4
		17	10
Porcelain	—	3	
	Porcelain Liner M	20	

TC : Thermal cycles

### Physical properties

Properties	Unit	Value
Flow	mm	55
Film thickness	μm	12
Compression strength	MPa	190
Flexural strength	MPa	57
Solubility in water	μg/mm <sup>3</sup>	12.4
Radiopacity	% Al	120
Easiness to remove excess	—	same as glass ionomers

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