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INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Protective Coat for Dentin

Hybrid Coat

FOR DENTIST USE ONLY

Hybrid Coat

1. What is Hybrid Coat?

Hybrid Coat is a single-bottle light-curing coating material for exposed dentin. It can be used to seal and protect the adhesive interface between tooth structure and resin composite or indirectly fabricated restorations. Hybrid Coat is also suitable for surface coating tooth surfaces that have suffered attrition and wear. It contains 4-META*, a high performance adhesive monomer that decalcifies tooth substrate and penetrates through the smear layer to form a hybrid layer. This layer forms a fortified resin coat on dentin, which blocks mechanical and biological stimuli to dentin-pulp-complex.

* 4-methacryloxyethyl trimellitate anhydride

2. Contents of the kit:

- ① Base ----- 5mL
- ② Coat Sponge ----- 1 (approx.160 pcs)
- ③ Plastic Dispensing Dishes ----- 5

Sold separately :

- Base ----- 5mL
- Coat Sponge ----- 1 (approx.160 pcs)
- Plastic Dispensing Dishes ----- 30

3. Precautions: Read all instructions thoroughly before use.

3.1 Safety

① Sensitivity

Hybrid Coat should not be used by clinicians or on patients who are methacrylate monomer-sensitive. If signs of irritation such as redness appear, stop using it immediately and consult a physician.

② Precautions

Hybrid Coat contains acidic methacrylate monomer. Avoid contact of the Hybrid Coat components with soft tissue, skin or eyes. Avoid ingestion of the components. To avoid contact with mucosa or gingiva, use of rubber dam or a protector such as vaseline is recommended. Contaminated skin or mucosa should be wiped off immediately with alcohol and then rinsed thoroughly with copious amounts of running water. In some cases contaminated mucosa or gingiva may whiten, which will usually disappear in a few days. To minimize contamination during blowing or coating, place an evacuator near the treated tooth.

If signs of skin sensitization such as rash, redness, ulcer, swelling, itch, numbness on mucosa, gingiva or lips appear, stop using Hybrid Coat immediately and consult a physician.

If the Coat Sponge contacts wet skin, the odor may remain.

If Hybrid Coat enters the eye, immediately rinse the eye thoroughly with running water. The patient should then be examined by an ophthalmologist.

Clinicians should use dental protective gloves. Care should be taken to avoid contact with patients with the gloves which previously contacted Hybrid Coat.

Base contains approx. 40% acetone. Inhalation of a large amount of concentrated vapor may cause headache: therefore the area should be well ventilated. In case of such severe inhalation, move to fresh air.

③ Pulp protection

If the preparation approaches the pulp, apply a protective liner or base.

④ Flammability

Base is flammable. **AVOID OPEN FLAME** in use and storage.

⑤ Applications

Use Hybrid Coat only for the applications recommended in the instructions.

3.2 To get the best results

① Coat sponge

Coat Sponge contains a bonding promoter and is an important component of Hybrid Coat. Always apply Hybrid Coat using Coat Sponge and **do NOT substitute any other sponges.**

② Storage

Store Hybrid Coat in a dark location between 4°C and 30°C, preferably in a refrigerator, but **do NOT freeze Hybrid Coat.** If frozen, phase separation will occur. High temperature, direct sunlight and high humidity should be avoided. Coat Sponge is slightly hydroscopic.

③ Volatility

The Base is highly volatile. **Firmly recap the bottle immediately after use.**

④ Contamination

To avoid contamination, do not mix the components with those of other materials, nor mix the bottle caps. Particularly, do not mix Coat Sponge with any other sponges. After use, dispose of leftover material. Once dispensed, do not return the materials to their original container.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

4. How to use Hybrid Coat

1) Coating indirect restoration

To avoid contact with mucosa or gingiva, use of rubber dam or a protector such as vaseline is recommended.

① Surface preparation

Prepare tooth using an instrument such as fine-grit point to minimize smear layer formation. After tooth preparation, wash and blow-dry the surface.

Note : The bond strength to enamel can be further improved by pretreatment of enamel for 30 seconds with a phosphoric etchant such as Super-Bond C&B Red Activator but do NOT use silica-thickened phosphoric etchant which will compromise the bond strength. Fluorinated enamel should be always etched as above.

② Application of Base

a) Dispensation

Dispense one drop (up to 3 drops for a larger surface) of Base in the well of a Plastic Dispensing Dish and firmly recap the bottle immediately.

Stir the expressed Base in the well for a few seconds with the Coat Sponge supplied in the Hybrid Coat kit.

Note : The mixture in the mixing well should be used within 3 minutes to minimize evaporation.

b) Application

After stirring, apply the mixture onto the prepared tooth surfaces and keep them moist for 10 - 20 seconds. To avoid the excessive mixture contact with gingiva, use of a marketed brush applicator is strongly recommended.

Air-blow the solvents for 5 - 10 seconds. Apply a gentle stream of air for 5 - 10 seconds where the mixture lays thick prior to the strong air for 5 - 10 seconds.

Note : Using vacuum suction is recommended. Place the vacuum tip close to cover the treated surface so that the Hybrid Coat does not splash intraorally during the air-blow.

③ Light-curing

Irradiate with a visible light curing unit.

Light	Light Strength	Irradiation Time
Halogen	120mW/cm ² or above	ca. 5 - 10 sec
Plasma	2,000W/cm ²	ca. 3 sec
LED	600 - 800mW/cm ²	ca. 5 sec

Note : Hybrid Coat is compatible with all currently available curing lights. With a fast light, the irradiation time can be shortened according to the manufacturer's instructions.

④ Restoration

a) Restoration with a resin composite

Place a resin composite on the coated surface and light-cure according to the manufacturer's instructions.

b) Restoration with an indirectly fabricated prosthesis

To remove the air-inhibited layer, make sure to wipe the surface of the light-cured Hybrid Coat with alcohol on a cotton pellet. Where necessary, re-shape the margin.

Take an impression using an alginate/agar impression material.

Note : Make sure if the polymerization is adequate in case a silicone impression materials is used.

Temporary material or temporary crown should be used to protect the coated surface. If a temporary resin is used, apply a separating agent such as Super-Bond SEP beforehand to avoid inadvertent bonding to the tooth.

Note : Do not use resin-based or eugenol-based temporary materials since they may compromise the polymerization of luting cements.

Remove the temporary, wash and dry the cavity/prepared surface.

Apply a luting agent of your choice and place the prosthesis. Resin cement such as Super-Bond C&B is recommended.

2) Surface Coating of Attrition and Wear

To avoid contact with mucosa or gingiva, use of rubber dam or a protector such as vaseline is recommended.

① Surface preparation

Clean the tooth surfaces to be coated. Rinse thoroughly and dry.

Note : Non-cut enamel should be etched for 30 seconds with a phosphoric etchant such as Super-Bond C&B Red Activator. Do NOT use silica-thickened phosphoric etchant which will compromise the bond strength. Fluorinated enamel should be always etched as above.

② Application of Base

Place a retraction cord around the gum line to eliminate Hybrid Coat from entering the gingival sulcus.

a) Dispensation

Dispense one drop (up to 3 drops for a larger surface) of Base in the well of a Plastic Dispensing Dish and firmly recap the bottle immediately.

Stir the expressed Base in the well with a piece of Coat Sponge for a few seconds.

Note : The mixture in the mixing well should be used within 3 minutes to minimize evaporation.

b) Application

After stirring, apply the mixture onto the prepared tooth surfaces and keep them moist for 10 - 20 seconds. To avoid the excessive mixture contact with gingiva, use of a marketed brush applicator is strongly recommended.

Air-blow the solvents for 5 - 10 seconds. Apply a gentle stream of air for 5 - 10 seconds where the mixture lays thick prior to the strong air for 5 - 10 seconds.

Note : Using vacuum suction is recommended. Place the vacuum tip close to cover the treated surface so that the Hybrid Coat does not splash intraorally during the air-blow.

③ Light-curing

Irradiate with a visible light curing unit. Use the irradiation time listed in the table in "③ Light-curing" of "1) Coating indirect restoration".

Note : Should the attrition or wear be severe, restoration using resin composite may follow.

Note : Should the treatment procedure end with the surface coating with Hybrid Coat, the coated surface should be wiped off with a cotton pellet sucked in alcohol to remove air-inhibited layer. The excessive alcohol should be squeezed out of the cotton pellet before use.

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*4-methacryloxyethyl trimellitate anhydride

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- ①Base ----- 5mL
 - ②Coat Sponge ----- 1통 (약 160 개)
 - ③Plastic Dispensing Dishes ----- 5 개
- [별매품]**
- Base ----- 5mL
 - Coat Sponge ----- 1통 (약 160 개)
 - Plastic Dispensing Dishes ----- 30 개

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