EN **FR** **DE** **IT** **ES** **TR**

EN

INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Flowable Radiopaque Hybrid Light-curing
Resin-based Dental Restorative Material

Metafil Flo Medium-Flow

FOR DENTIST USE ONLY

1. What is Metafil Flo Medium-Flow?

Metafil Flo Medium-Flow is a highly-flowable, visible-light-cure, tooth-shade, radio-opaque composite resin based on monomer combination and containing inorganic fillers and an organic filler which is reactive with the matrix. The inorganic filler loading is 42.4vol% with a particle size range of 0.01 to 1µm.

Indications:

Metafil Flo Medium-Flow is used for a restoration of defective teeth in the oral cavity. Do not use for a restoration of root canals.

2. Types

- Metafil Flo Medium-Flow ----- 1.5mL
Shade: A1, A2, A3, A3.5, B2, CV, D-BROWN, OPAQUE, U-OPAQUE
- 20G Needle Tip ----- 10 (with 1 Needle Cap)

3. Components

Methacrylate monomers, Organic filler, Inorganic filler, Aromatic amine, etc.

4. Precautions

- This product contains substances that may cause allergic reactions.
- Metafil Flo Medium-Flow should NOT be used either by clinicians or on patients who are methacrylate-monomer sensitive.
- Avoid contact of Metafil Flo Medium-Flow with soft tissues, skin or eyes. Contaminated skin or mucosa with Metafil Flo Medium-Flow should be wiped off immediately with alcohol and then rinsed thoroughly with running water. If Metafil Flo Medium-Flow enters an eye, rinse the eye with running water immediately. The patients should then be examined by an ophthalmologist.
- Avoid accidental ingestion while applying and rinsing the material. If ingested, the patient must seek medical attention.
- Use a dental rubber dam to isolate the treatment area from contamination.
- Needle Tip is for one-time use only. It is not autoclavable or chemically sterilizable. Dispose it after each patient to avoid cross-contamination.
- Use the Needle Tips only for this material.
- Dentist should wear rubber or PVC dental gloves.
- In case dispensed material is not used immediately, cover it to protect the material from ambient light to prevent polymerization.
- Use Metafil Flo Medium-Flow only for the applications recommended in this leaflet.
- If necessary, the pulp must be indirectly capped locally for deep cavities with only a minimal layer of dentin.
- Product lifetime can be shortened depending on the circumstances of use. Read instructions of relevant products carefully before use.

5. Storage

Store Metafil Flo Medium-Flow at a room temperature between 1-30°C.
Avoid high temperatures above 30°C, high humidity and direct sunlight. Refrigeration is not necessary. However, if stored in a refrigerator, then allow the material to warm to room temperature prior to use.
Shelf life: 3 years in a room temperature. See outer package for the expiration date. Do not use after expiration date.

6. Disposal

Dispose of empty containers, etc. in accordance with local regulations.

7. How to use Metafil Flo Medium-Flow

Adhesive Bonding System :

Apply light-cured bonding system according to manufacturer's instructions. We recommend the use of Hybrid Bond from Sun Medical Co., Ltd.

Placement and Curing :

- Remove the cap by turning it counter-clockwise.
- Screw the Needle Tip tightly to the syringe.
- Place dispensed material into prepared cavity previously treated with the bonding agent.
- Light-cure using dental visible light-curing device according to the manufacturer's instructions.

Type	Curing time (sec)	Depth of Cure (mm)	
		A1, A2, A3, A3.5, B2, CV	D-BROWN, OPAQUE, U-OPAQUE
Conventional halogen	20	1.5	1.0
	10	1.5	1.0
Fast halogen	20	2.0	1.0
	3	1.5	1.0
Plasma arc	10	1.5	1.0
	20	2.0	1.0

- For best handling, use Metafil Flo Medium-Flow at room temperature.
- If the cavity preparation is large and/or deep, incremental filling techniques should be employed.

Polishing :

Metafil Flo Medium-Flow will polish quickly and easily to a high luster. Use of most discs or polishing paste systems are acceptable.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

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