





Metafil Flo is a highly-adaptable resin composite. It contains inorganic fillers and a unique organic filler which is very reactive with the matrix resin. This assures superior handling and highly wear-resistant restoration. Metafil Flo is radio-opaque, with low shrinkage and deep depth of cure, it is suitable to use for base build up material as well.

Typical Applications

Restoration of Attrition









Metafil Flo has high wear resistance for restoring attrition.

Cervical Lesion









Higher flexural strength of Metafil Flo will decrease the dislodging of restoration.

Restoration of Proximal Cavity









As a guick polish. Metafil Flo is highly esthetic and ideal for anterior restoration.

Restoration of Sallow Enamel Lesions / Sealants









Metafi Flo is highly wear-resistant and can be used for small enamel lesions / sealants.

Shades



















D-Brown Opaque U-Opaque