



Light-cure Orthodontic Adhesive

Orthomite LC

Easy to Use, Time-saving
Reliable Bonding Performance
Applicable to Any Brackets






Orthomite LC is a light-cure adhesive resin for orthodontic bonding of metal, ceramic and plastic brackets to tooth surface.


Easy to Use, Time-Saving

Through the innovative design of the base monomer, Orthomite LC Paste can bond to etched enamel well without priming.

Reliable Bonding Performance



Clinical Procedure



Remove excess adhesive before curing.

Metal brackets : both interproximal sides 10 sec. Ceramic brackets and plastic brackets : 10 sec. through brackets.